1

(227) 4-AMINOPHENOL IN ACETAMINOPHEN-CONTAINING DRUG **PRODUCTS**

INTRODUCTION

This chapter provides a procedure and acceptance criterion (limit) to control the principal degradation product of acetaminophen, 4-aminophenol, an impurity that can form by hydrolysis of acetaminophen.

SOLUTION PREPARATIONS

All solution preparations that contain acetaminophen or 4-aminophenol should be protected from light and should be stored only for as long as can be supported by solution stability data acquired during verification under actual conditions

Buffer: 4.0 g/L of sodium citrate dihydrate and 1.5 g/L of anhydrous citric acid, in water

Diluent: Acetonitrile and *Buffer* (10:90)

Solution A: 10 mM phosphate buffer prepared as follows. Add 0.60 g of monobasic potassium phosphate and 0.82 g of anhydrous dibasic sodium phosphate to a 1-L volumetric flask. Dissolve and dilute with water to volume to a pH of 7.0.

Solution B: Water Solution C: Acetonitrile Mobile phase: See Table 1.

Table 1

Time (min)	Solution A (%)	Solution B (%)	Solution C (%)
0	90	5	5
5	90	5	5
7	10	10	80
7.1	90	5	5
10	90	5	5

Standard stock solution: 25 µg/mL of USP 4-Aminophenol RS in Diluent. Prepare fresh in conjunction with the other solution preparations described below. Discard after 4 h or as supported by solution stability data. [NOTE—See Chromatographic Adjustments, solution stability, below.]

System suitability solution: 2.5 µg/mL of USP 4-Aminophenol RS in *Diluent*, from the *Standard stock solution*Sample stock solution: Nominally 10 mg/mL of acetaminophen from a suitable quantity of drug product in *Diluent*. [Note—Either component of the Diluent may be introduced to the drug product first, followed by addition of the other component to maintain the proportions of acetonitrile and Buffer and to achieve the appropriate final volume defined for the *Diluent*.]

[NOTE—It is recommended that the Sample solution and Standard solution be prepared concurrently within a narrow

window of time (e.g., 30 min) for each drug product sample.]

Standard solution: Add 25.0 mL of the Sample stock solution and 15.0 mL of the Standard stock solution to a 50-mL volumetric flask, and dilute with Diluent to volume. Pass through a suitable filter of 0.45-µm pore size, discarding the first 3 mL of filtrate.

Sample solution: Add 25.0 mL of the Sample stock solution to a 50-mL volumetric flask, and dilute with Diluent to volume. Pass through a suitable filter of 0.45-µm pore size, discarding the first 3 mL of filtrate.

CHROMATOGRAPHIC METHOD

Chromatographic system

(See Chromatography (621), System Suitability.)

Mode: LC

Detector: UV 300 nm

Column: 4.6-mm × 15-cm; 5-µm packing L85

Column temperature: 30° Flow rate: 1 mL/min Injection volume: 10 μL

System suitability

Samples: System suitability solution and Standard solution

[Note—The typical retention time for 4-aminophenol is about 4.2–5.3 min.]

Suitability requirements

Resolution: NLT 1.0 between 4-aminophenol and the nearest peak, Standard solution

Tailing factor: NMT 1.5 for the 4-aminophenol peak, Standard solution

Relative standard deviation: NMT 5.0%, Standard solution

Signal-to-noise ratio: NLT 20 for the 4-aminophenol peak, System suitability solution

Samples: Standard solution and Sample solution

Inject the Sample solution and Standard solution for each drug product sample sequentially, i.e., back-to-back.

Calculate the percentage of 4-aminophenol (C_6H_7NO) relative to acetaminophen in the portion of drug product taken:

Printed by: Nguyen Nhan

2

Official Date: Official as of 01-May-2018

Document Type: GENERAL CHAPTER

@2021 USPC

Result = $[r_{U}/(r_{S} - r_{U})] \times (W_{S}/W_{U}) \times 100$

= peak response of 4-aminophenol from the Sample solution r_{U} = peak response of 4-aminophenol from the Standard solution

 W_{s} = amount of USP 4-Aminophenol RS added to the Standard solution (mg)

= amount of acetaminophen in the Sample solution (mg) W_{U}

Acceptance criteria (unless otherwise stated in the monograph): NMT 0.15% of 4-aminophenol relative to acetaminophen

CHROMATOGRAPHIC ADJUSTMENTS

The retention time of 4-aminophenol can be tuned to achieve specificity for a given product matrix. This allowance supersedes provisions in $\langle 621 \rangle$ for adjusting chromatographic conditions and is intended to provide a measure of flexibility when needed. Suggestions for changing 4-aminophenol retention are given in Table 2. The use of a ternary mobile phase system affords ready changes to the ionic strength (water from Solution B) and organic strength (acetonitrile from Solution C), but this can be simplified to a binary mobile phase system.

Table 2

Condition Change	Change in 4-Aminophenol Retention	
Increase in organic strength (Solution C)	Decreases 4-aminophenol retention	
Decrease in pH (Solution A)	Increases 4-aminophenol retention	
Increase in ionic strength (Solution B)	Decreases 4-aminophenol retention	
Increase in column temperature	Increases 4-aminophenol retention	

Adjustments to the chromatographic procedure may require verification or validation. See Validation of Compendial Procedures (1225) and *Verification of Compendial Procedures* (1226) for guidance. Adjusted chromatographic conditions must meet all system suitability requirements.

Solution stability must be verified under actual conditions of use to ensure 4-aminophenol is stable in the Sample solution and the Standard solution as evidenced by NMT ±10% change of the 4-aminophenol peak areas.

• USP REFERENCE STANDARDS (11)

USP 4-Aminophenol RS

4-Amino-1-hydroxybenzene.

C₆H₇NO

109.13